

BUFFER SOLUTIONS

ST-MS Buffer Solution



FEATURES

- Stabilized liquid solution

FIELD APPLICATION:

The precision and reliability of Conductivity measurement is determined by the buffer solution used for calibrating the probe. The special double plug container ensures that a new unpolluted buffer is always available.

	Solution	Value	Quantity	code
ST MS 8	Conductivity	84 μ S/cm at 25°C	500 ml	9900122018
ST MS 14	Conductivity	1423 μ S/cm at 25°C	500 ml	9900122019
ST MS 128	conductivity	12880 μ S/cm at 25°C	500 ml	9900122020