

BUFFER SOLUTIONS

ST- pH Buffer Solution



Features:

- Stabilized liquid solution

Field Application:

The precision and reliability of pH measurement is determined by the buffer solution used for calibrating the probe. The special double plug container ensures that a new unpolluted buffer is always available.

	Solution	Value	Quantity
ST PH 4	pH	4,00 pH at 20°C	70 ml
ST PH 7	pH	7,00 pH at 20°C	70 ml
ST PH 10	pH	9,22 pH at 20°C	70 ml

Expire date indicated on the label product

Life time after opening is 30 days only, keep 20°C and not light sun directly

Life time after opening is 180 days only, keep 6°C and not light sun directly